



Puroxi



UltraSound Gel

Use for Ultrasound & Laser Machines

"Clear"



Hypoallergenic - Water Soluble - Alcohol Free

Keep Away From Sunlight , Not to be Swallowed

Produced Under GMP , Temperature: -10°C to $+50^{\circ}\text{C}$

Ingredients : Deionized Water , Carbomer , Triethanolamine , EDTA ,
Monopropylene Glycol , Preservative

Manufactured by SADARA GROUP

+961 71 9131 21

+961 8 9130 96



C.R.:4004720
Ind.R.:2535
Inv.R.:6153/T



500 ml